

Infusion Handles



Product Name: Infusion Handle w/ On / Off Control

Product Description: A flow on/ off control infusion handle for use with roller pump devices. Handle Diameter 15/16" (Stainless Steel Control) 3/4" Handle. Valve / Female Luer Lock / Male Luer Lock SS. Available Black or Blue.

Black Part #: A005

Blue Part #: A006




Product Name: Infusion Handle with on / off stopcock control

Product Description: Infusion Handle w / stopcock control, Handle Diameter 3/4" (Small) 1" (Large), Valve: Brass w Electroless Nickel Plating. Female and Male Luer Lock SS.

Small Part #: A066

Large Part #: A065



Product Name: Infusion Handle (Small & Large)

Product Description: Standard Infusion Handle, No Control, ideal for infusion needle. Handle Diameter- Sm (3/4") Lg (1"), Female and Male Luer Lock SS.

Small Part #: A004

Large Part #: A003



Product Name: Infusion Handle with Swivel Luer Nut (Small & Large).

Product Description: Orientation Handle, ideal for positioning the cannula port to handle thumb rest. Handle Diameter- Sm (3/4") Lg (1") Female and male swivel luer lock SS.

Small Part #: A064

Large Part #: A063



Product Name: Luer Lock Adjustable Handle (Small and Large)

Product Description: Orientation handle ideal for positioning cannula port to handle thumb rest. Lock by tightening the cap. Handle Diameter 1" (Large), female and male luer Lock SS.


Large Part #: A062



Product Name: Luer Lock Adjustable Aspiration Handle (Small, Large)

Product Description: Luer Lock adjustable handle. Ideal for rotating and locking position desired of the cannula port to thumb rest. Handle Diameter 3/4" (small) 1" (large), Male Luer Lock SS.

Large Part #: A061




Product Name: Aspiration Handle (Small, Large)

Product Description: Standard Aspiration Handle, Handle Diameter 3/4" (small), 1" (large), male luer lock SS.

Small Part #: A007

Large Part #: A008



Product Name: Aspiration Handle w/ Double Luer Lock 1" (large)

Product Description: Standard Aspiration with Double Luer Lock, Handle Diameter 3/4" (small), 1" (large). Male and female luer lock.

Large Part #: A009